

Belscan Continental BV Stokkelaar 4 9160 Lokeren

Belgium

Tel: +32 9 340 56 26

Email: info@belscan.be

					Home: www.l	Belscan.be
PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9001325	9000894	9000895	9000896	9000897	
Description	Nitrile Gloves		HY@PRO	N-3,2	AQL 1,5	
	Nitrile gloves	blue powder	free M3.2 10	Opcs/dispense	er box HY@PRO	

### PRODUCT PACKAGING IMAGE



PRODUCT DESCRIPTION										
Material			Nitrile							
Colour		Blue								
Characteristics		Powderfree		Non Sterile	Ambidextrous					
Surface		Finger Textured								

SIZES										
	XS (5-6)	S (6-7)	M (7-8)	L (8-9)	XL (9-10)	XXL (10-11)	XXXL (11-12)			
Width		80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	115 ± 10 mm					
length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm					

THICKNESS									
	Cuff	Palm	Fingertips	Overal weight (M)					
Thickness	0,08mm	0,10 mm	0,10mm	3,2 gram					

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9001325	9000894	9000895	9000896	9000897	
Description	Nitrile Gloves					
	Nitrile gloves 100pcs/dispe	•		4		The state of the s


STANDARDI:	SATION EN388	3						
EN 388 Mechanical risks	Abrasion resistance	Blade cut resistance	Tear resistance	Punture resistance	Blad cut resistance TDM-test	Impact test		EN 388:2016
Level								
STANDARDI	SATION EN374							
EN 374-1	The gloves must not	change their shape to	oo much, nor must th	ey break or stick.			<b>~</b>	
EN 374-2	The material must no	ot contain any perfora	ations				<b>~</b>	
EN 374-3	Chemicals must not	nemicals must not pass through the glove material						
EN 374-4	Tests are carried out on the duration and the amount of chemical substances that penetrate the glove.							
EN 374-4	Chemical		Code letter	Permeation level		Degradation		
	Sodium Hydroxide		К	6		-33,50%	<b>~</b>	
	Hydrogen Peroxide		Р	4		22,10%	<b>~</b>	
	Formaldehyde		Т	6		-0,003	<b>~</b>	
							$\circ$	
							$\circ$	
Accelerators							$\circ$	
EN 374-5	The glove is tight aga	ainst microorganisms	(viral, bacteria and fu	ıngi).				EN ISO 374-5-2016
Microorganism tightness	Testing method		~	VIRUS				
STANDARDI	SATION EN165	523						
EN16523-1		aterial resistance to pe	•	als - Part 1: Permeatio	on by potentially haza	rdous liquid	<b>V</b>	

STANDARDISATION EN455									
EN 455 medical gloves for single use	The glove meets the requirements according to EN 455-1, EN 455-2, EN 455-3, EN 455-4	<b>~</b>	EN455						
EN 455-1	The glove has an AQL $< 1.5$ in regards to the water retention test (sampling inspection in acc. to ISO 2859-1, general Inspection Level I)	<b>~</b>	EN455-1						
EN 455-2:	The tensile strength, thickness and fit is tested with a minimum standard of 6 Newtons.	<b>&gt;</b>	EN455-2						
EN 455-3	This biological evaluation tests for the presence of harmful substances such as powder and micro-organisms.	<b>✓</b>	EN455-3						
EN 455-4	The product tested shelf life	<b>~</b>	EN455-4						
STANDARDISATION EN21420									
ISO EN 21420 Protective gloves	The glove meets the requirements according to EN 420 / ISO EN 21420	<b>~</b>							

Rev-nr: 001 Date: 24/06/2024

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9001325	9000894	9000895	9000896	9000897	
Description	Nitrile Gloves					
	Nitrile gloves blue powder free M3.2 100pcs/dispenser box HY@PRO					M.

REGALUTOR'	Y AFFAIRS								
PPE-Regulation (EU) 2016/425			Category III						
MD-Regulation (EU) 2017/745									
CE Regulation		SATRA 2777							
EUROPEAN REPRESENTATIVE							$\circ$	EC REP	
Food Contact (EG) 1935/2004	acidic foods	aqueous foods	fatty foots	alcoholic foods	dry foods		>	Ω"	
LATEX	Does not contain late	Does not contain latex							
Storage conditions	Store in original packaging in a dry and dark place at 10°C to 30°C. Refer to guidelines of storage of rubber product as decribed in ISO2230:2002. Ensure that storage is kept cool, dry and dust-free, avoid ventilation and storage close to phocopy equipment. Copper-ions discolour the glove.  Shelf life (months) 36 Protect gloves against ultraviolet light sources, as sunlight and oxidizing agents. Storage above 30°C will lead to accelerated aging and should be avoided.							***	
Single/double/Reus able	Examination glove fo	or single use					<b>~</b>	2	
Sterile / Non Sterile	Non Sterile glove						>	NOR	
ASTM F739	Symbol used for SW	chemotherany tested	gloves that have bee	n tested either hy AST	TM F739 or by ASTM	D6978		TESTED FOR USE WITH	
ASTM D6978	Symbol used for Sw	enemotherapy tested	I gloves that have bee	in tested entirel by As	TIVIT 735 OF BY ASTIVIT	50370	)	CHEMOTHERAPY DRUGS	
REACH Compliant	REACH is a regulation of the European Union, adopted to improve the protection of human health and the environment from the risks that can be posed by chemicals, while enhancing the competitiveness of the EU chemicals industry.							( jestvij)	
CE	CE conformoty							(€	
Intrastat number	4015190000								
Made in	CN								



## LOGISTIC DATA SALES UNIT

General information		Material	Carton					
Size	pieces	EAN Code	REF code	Length (mm)	Width (mm)	Heigth (mm)		Gross weight
XS	100	5412276278659	9001325	220	120	60	290	340,9
S	100	5412276201947	9000894	220	120	60	290	340,9
М	100	5412276202043	9000895	220	120	60	320	370,9
L	100	5412276202111	9000896	220	120	60	380	430,9
XL	100	5412276202166	9000897	220	120	60	450	500,9
XXL							·	·

## LOGISTIC DATA OUTER PACKAGING

General information		Material	Carton					
Size	pieces	EAN Code	REF code	Length (mm)	Width (mm)	Heigth (mm)		Gross weight
XS	10	5412276278666	9001325	320	250	230	3,11	3,41
S	10	5412276241417	9000894	320	250	230	3,409	3,7054
М	10	5412276241431	9000895	320	250	230	3,709	4,0054
L	10	5412276241455	9000896	320	250	230	4,309	4,6054
XL	10	5412276241691	9000897	320	250	230	5,009	5,3054
XXL								

# LOGISTIC DATA PALLET

General information		Pallet: EURO-PALLET						
Size	Cartons per layer	Layers per pallet	Netto height	Gross Heigth (incl. Euro pallet 14,4cm)	Gross weight (incl. Euro pallet 25kg)			
XS	11	6	23	152,4	250,06kg			
S	11	6	23	152,4	269,56kg			
М	11	6	23	152,4	289,36kg			
L	11	6	23	152,4	328,96kg			
XL	11	6	23	152,4	375,16kg			
XXL								

Rev-nr: 001 Date: 24/06/2024

PRODUCT DATA SHEET						
	XS	S	М	L	XL	XXL
Our reference	9001325	9000894	9000895	9000896	9000897	
Description	Nitrile Gloves					
	Nitrile gloves blue powder free M3.2 100pcs/dispenser box HY@PRO					

#### **Declaration of conformity Food**



Food Compliance Declaration - Belscan

Date: 04-07-2022 Owner: Q-Team Code: I.02.09.03.DOC.10 Version: 001 Authorizer: QA

# Declaration of compliance for materials and articles intended to come into contact with food

The Supplier, or his authorized representative who is established in the community:

Name: Belscan Continental BV

Address: Stokkelaar 4-9160 Lokeren, Belgium

Declares that the following product(s):

Ref. 9001325-9000894-95-96-97: Nitrile gloves blue powder free M3,2 100pcs/dispenser box HY@PRO

is in accordance with Regulations (EC) No 2023/2006 (Good Manufacturing Practices)<sup>1</sup>, and No. 1935/2004 (Non-Plastic Materials)<sup>2</sup>.

For use in the following circumstances

- type of food stuff intended for coming into contact with materials / objects: All types of food
- Storage time and temperature of the material/object: 3 years- between 5° c and 30 °c
- Possible treatment for the material/object: not suitable for heat treatment
- Ratio surface/volume: 10 mg/d²

(Appropriate information on all substances for which there are restrictions, at the level of the EU as well as at the Belgian level, so that all future users can comply with those restrictions. In the absence of any national or European regulation, all information on international restrictions, norms or guide values has to be provided (Council of Europe, WHO, Codex Alimentarius...))

Completed at,

Lokeren

Radea Ahaddour Quality department For Belscan Continental BV Stokkelaar Gokfurn Bet Undows 867

Date: 24/06/24

Rev-nr: 001 Date: 24/06/2024